

AFLP Protocol

28 October 2011

Anthony J Geneva

Modified from:

Vos, P., R. Hogers, M. Bleeker, M. Reijans, T. Lee, M. Hornes, et al. 1995. AFLP: a new technique for DNA fingerprinting. *Nucleic Acids Research* 23:4407-4414.

and

<https://www.msu.edu/user/hazensam/aflp/AFLPprotocolMSU.html>

Oligos

Adaptors

EcoR1adapterA CTCGTAGACTGCGTACC
EcoR1adapterB AATTGGTACGCAGTCTAC
Mse1adapterA GACGATGAGTCCTGAG
Mse1adapterB TACTCAGGACTCAT

Pre Selective Primers

Prese1E GACTGCGTACCAATTCA
Prese1M GATGAGTCCTGAGTAAC

Mse Selective Unlabeled Primer

SelampM47 GATGAGTCCTGAGTAACAA
SelampM51 GATGAGTCCTGAGTAACCA
SelampM55 GATGAGTCCTGAGTAACGA
SelampM59 GATGAGTCCTGAGTAACTA
SelampM49 GATGAGTCCTGAGTAACAG
SelampM53 GATGAGTCCTGAGTAACCG
SelampM61 GATGAGTCCTGAGTAACTG

EcoR1 Selective 5'-Labeled Primer

SelAmpE1 GACTGCGTACCAATTCAGG
SelAmpE2 GACTGCGTACCAATTCATC
SelAmpE3 GACTGCGTACCAATTCACA

Dilute all stocks to 100 μ M and store -20°C.

Adaptor Construction

Component	Amount (ml)
Mse1adapterA (100 μ M)	7.5
Mse1adapterB (100 μ M)	7.5
dH ₂ O	135
Total (5μM)	150

Component	Amount (ml)
EcoR1adapterA (100 μ M)	75
EcoR1adapterB (100 μ M)	75
Total (50μM)	150

Cycler Protocol

95°C 5 mins

70X 94°C 30 secs/ -1°C per cycle

4°C hold

Store at -20C

Digestion

Component	Amount (μ l)	Per plate	/2
dH ₂ O	13.5	1417.5	708.75
10X NEB buffer #4	2	210	105
NEB EcoR1 (20 units/ μ l)	0.25	26.25	13.125
NEB Mse1 (10 units/ μ l)	0.05	5.25	2.625
NEB 100X BSA	0.2	21	10.5
DNA (50ng/ μ l)	4		
Total	20		

Cycler Protocol

37°C 180 mins

60°C 15 mins

Ligation

Component	Amount (μ l)	Per plate
dH ₂ O	5.8	609
10X T4 ligase buffer (with ATP)	2	210
EcoR1 adapter (5 μ M)	1	105
Mse1 adapter (50 μ M)	1	105
T4 DNA ligase (400 units/ μ l)	0.2	21
Entire Digestion Product	20	
Total	30	

Incubate overnight at 37°C.

Run 4 μ L of Ligation product on a 1.5% agarose gel. A diffuse smear (or sometimes distinct bands) should be visible between ~200-1000 bp. If these appear, perform a full of plate preselective amplifications. If not, test a small number of samples via preselective amplification and continue if these are successful.

Preselective Amplification

Component	Amount (μ l)	Per plate
dH ₂ O	22.8	2394
E Primer (10 μ M)	5	525
M Primer (10 μ M)	5	525
MgSO ₄	5	525
10x BioBasic Buffer	5	525
dNTPs	5	525
Taq	0.25	26.25
Ligation Product	2	
Total	50.05	

Cycler Protocol

94°C 120 secs

94°C 60 secs

26x 56°C 60 secs

72°C 60 secs

72°C 60 secs

4°C hold

Run 4 μ L of Preselective Product on a 1.5% agarose gel. A DNA smear or distinct bands should appear between 50-500 bp. Store product at -20°C.

Selective Amplification

Component	Amount (μ l)	Per Plate
H ₂ O	11.4	1197
Buffer	2.5	262.5
dNTPs	2.5	262.5
Mse1 selective primer (2 μ M)	2.5	262.5
EcoR1 labeled selective primer (2 μ M)	2.5	262.5
MGSO ₄	2.5	262.5
Taq	0.125	13.125
Template	1	
Total	25.025	

Cycler Protocol

94°C 120 secs

94°C 30 secs

12x 65°C 30 secs/ -1°C per cycle

72°C 60 secs

94°C 30 secs

23x 56°C 30 secs

72°C 60 secs

72°C 60 secs

4°C hold

Run 4 μ L of Selective Product. Bands should appear between 50-500 bp. Store product at 4°C in the dark by wrapping in foil.

Analysis

TBA